

Fast characterisation and detection platform for ultra-thin monolayer semiconductors

Glen Pearce and Lachlan Catto

Supervisor: Dr. Matt Broome, CSIRO


Academic Mentor: Dr. Daniel Terno, MQ

1. What is 2D material?
2. Software for automatic detection
3. User-friendly
4. Benefits of automation




Motivation

- 3D Semiconductors are a well researched technology
- Narrow this down to 2D devices - monolayers
- Efficiency!



<https://www.sciencedirect.com/science/article/abs/pii/S0378517318307002>




mpohodzhay / shutterstock.com




<https://www.linkedin.com/pulse/how-focus-business-efficiency-productivity-year-robert-ford/>

Background

- Transition-metal Dichalcogenide (TMD) Semiconductors
- Many TMD monolayers have direct (optically excitable) band gap



- Good candidates for transistors, and optical devices



Thin layers are obtained by:
Mechanical Exfoliation (ME)



Methods – Overall Picture

Fluorescence Microscopy to detect only the monolayers, as since they have a direct band gap, thus can be excited optically.



Motorised stage and scientific camera to move and image the sample.

Post-processing with OpenCV to detect regions of fluorescence.

Analysing images to determine location, area, quality.

Displaying these results in an accessible format.




The monolayers of interest could then be marked for further study.



Methods – Overall Scan and Detection Logic



Serpentine search over the whole sample:

1. Move the microscope stage a precalculated distance, just enough to have some overlap in the image
2. Take an image, and place it onto a larger canvas, based on where the image was taken and blending the edges



Methods - Overall Scan and Detection Logic

1. Downscale & blur to reduce noise and speed up processing
2. Convert to grayscale, with a bias towards the fluorescence colour.
3. Apply a threshold to remove all low intensity pixels
4. OpenCV - *Find Contours* - to outline the bright areas (hopefully monolayers)
5. Use the contour data to approximate area, location and to save individual Monolayer images for analysis



Methods – Issues we had to overcome

- Camera Resolution
 - Memory Usage
 - CPU usage
- Camera rotation & view rotation
- Exposure
- GUI responsiveness
- Processing speed
- Thorlabs Driver
- Threading issues
- Closing the hardware connections safely
- Focus (particularly in large scans)



Methods – Autofocus Problem

Measuring Focus


To optimise the focus, first we need to quantify it.

In our case, our parameter needs to be:

- Robust against noise
- Fast enough to calculate
- Exposure independent

Our final method utilised the Power Spectral Density of the image.

Well-defined edges need more high frequency components in the image.



The more high frequency components, the better the edge definition



Image
↓
2D Fourier Transform
(spacial freq. & amplitude)

↓
Radial Averaging
(counteracts noise and makes 1D)

↓
Power Spectral Density
(power vs spacial freq.)

↓
Linear fit (log-log plot)


↓
Gradient
(shallower gradient → more high freq. → better focus)




Methods – Autofocus Problem

Finding the focus

Now that the focus can be quantified, finding the optimal focus should be relatively easy.....




Automatic Monolayer Finder




Pos X 6.827 mm
Pos Y 10.990 mm
Focus Z 1.321 mm

Methods - Designing good software


Bad GUI



Good GUI



Threading




https://www.cs.uic.edu/~jbell/CourseNotes/OperatingSystems/4_Threads.html


- Easy to use GUI
- Effective Threading

Methods - Designing the GUI


Main Control Tab



Calibration Tab



Results Tab



Results Video

Automat

Automatic Monolayer Finder

Main Control

Calibration

Results & Analysis

50um

50um

Pos X

Pos Y

Focus Z

Pos X

Pos Y

Focus Z

6.827 mm

17.342 mm

1.300 mm

Top Left

Bottom Right

Begin Search

Top Left: ['1.58e+07', '1.70e+07']

Bottom Right: ['6.83e+06', '1.10e+07']

Scan Complete!

Move Controls

Top Left

Up

Top Right

Left

Centre

Right

Bottom Left

Down

Bottom Right

Focus +

Focus -

Enter Positions (mm):

X:





Y:

Z:

es

Results Images

Examples of stitched images and monolayers



Results

- Produced a reliable software
- Incorporated useful features
 - Auto Focusing
 - Auto Exposure
 - Progress bar
- Increased throughput
- Software produces a stitched image of the scan area
- Excellent detection (with some false positives ~5-20% depending on the contamination)
- Cat button

Approximately 10x

Condensing a day of manual work into about 30-60 mins




Easily removed when inspecting the final result

Conclusions

Takeaways:

- 2D materials are a single or few atoms thick
- Monolayer TMDs are a new and promising area of 2D material science
- Our software makes monolayer detection after mechanical exfoliation much easier



WS₂, an example of a monolayer TMD.
https://en.wikipedia.org/wiki/Tungsten_disulfide